



# FINA/ARENA Swimming World Cup 2009

DURBAN MOSCOW STOCKHOLM BERLIN SINGAPORE

21-22 November 2009

Event 23

Women's 400m Individual Medley  
400m 4-Nages Femmes

Final  
Finale



## Results

Résultats

	Record	Splits			Name	NOC Code	Location	Date
<b>WR</b>	4:25.06	29.60	1:02.28	1:36.63	2:10.12	BELMONTE GARCIA Mireia	ESP Rijeka (CRO)	14 DEC 2008
		2:46.86	3:24.51	3:55.13				
<b>WC</b>	4:26.40	29.29	1:01.73	1:36.22	2:09.84	BELMONTE Mireia	ESP Stockholm (SWE)	11 NOV 2009
		2:47.00	3:24.42	3:56.19				

Final

Event No. 23

Rank	Lane	Name	NOC Code	R.T.	Time	FINA Points
<b>1</b>	<b>4</b>	<b>MEAKLIM Kathryn</b>	<b>RSA</b>	<b>0.85</b>	<b>4:22.88</b>	1049 WR
50m (1) 28.84	100m (1) 1:01.68 32.84	150m (1) 1:36.73 35.05	200m (1) 2:10.57 33.84	250m (1) 2:46.62 36.05	300m (1) 3:22.75 36.13	350m (1) 3:53.48 30.73
<b>2</b>	<b>5</b>	<b>KARPEEVA Svetlana</b>	<b>RUS</b>	<b>0.91</b>	<b>4:38.35</b>	884
50m (2) 30.27	100m (2) 1:04.40 34.13	150m (2) 1:40.62 36.22	200m (2) 2:16.42 35.80	250m (2) 2:54.37 37.95	300m (2) 3:33.26 38.89	350m (2) 4:06.36 33.10
<b>3</b>	<b>3</b>	<b>KOH Ting Ting</b>	<b>SIN</b>	<b>0.83</b>	<b>4:54.81</b>	744
50m (3) 30.73	100m (3) 1:05.54 34.81	150m (3) 1:44.02 38.48	200m (3) 2:21.60 37.58	250m (3) 3:04.16 42.56	300m (3) 3:46.58 42.42	350m (3) 4:21.21 34.63
<b>4</b>	<b>6</b>	<b>LIM Meagan</b>	<b>SIN</b>	<b>0.78</b>	<b>5:03.48</b>	682
50m (4) 31.37	100m (4) 1:07.21 35.84	150m (4) 1:45.59 38.38	200m (4) 2:23.83 38.24	250m (4) 3:07.98 44.15	300m (4) 3:52.16 44.18	350m (4) 4:28.33 36.17
<b>5</b>	<b>7</b>	<b>CHAN Shi Min Marina</b>	<b>SIN</b>	<b>0.73</b>	<b>5:13.97</b>	616
50m (5) 31.59	100m (5) 1:09.14 37.55	150m (6) 1:50.10 40.96	200m (6) 2:30.07 39.97	250m (5) 3:16.23 46.16	300m (5) 4:02.50 46.27	350m (5) 4:38.67 36.17
<b>6</b>	<b>2</b>	<b>TAN Cheryl</b>	<b>SIN</b>	<b>0.80</b>	<b>5:16.77</b>	600
50m (6) 32.00	100m (6) 1:09.30 37.30	150m (5) 1:49.69 40.39	200m (5) 2:30.01 40.32	250m (6) 3:16.50 46.49	300m (6) 4:03.45 46.95	350m (6) 4:40.89 37.44

**Legend:**

R.T. Reaction time

WR World record

Timing & Data Handling by OMEGA

